

"Zephyr" - the ideal Transit Wheelchair.

Designed to take all the effort out of pushing wheelchairs, the "Zephyr" provides the perfect solution to wheelchair transfers.

With its simple, intuitive driving controls including a twist-grip throttle and direction control, variable speed governor and light, handlebar steering, the "Zephyr" provides a comfortable and safe transfer for the patient - and total safety for the operator.



The comfortable, high-back, office style chair is covered with easy-clean leather-look material with comfortable rounded armrests, while the flip-up footrests are fitted with removeable heel guards. A lap type seat belt is provided for added patient safety.

Safe, very stable, quiet and easy-to-use, the "Zephyr" will travel approximately 20kms between battery charges, allowing almost continuous use and making light work of heavy patient transfers over long distances in complete safety. Where Bariatric patients are transported, a heavy-duty version is available. The "Zephyr HD" is fitted with a 20" seat and is designed to carry patients up to 225kgs in complete safety.

Technical Specification:

Power System:	2 x 12 Volt, 24A/h Deep Cycle Lead Acid Gel batteries
Motor:	24 Volt, Heavy-duty electric motor and Gearbox
Drive system:	Heavy-duty, minimum-stretch chain and sprocket
Controller:	Dynamic High Capacity Programmable Controller
Drive wheels:	260 mm pneumatic, non-marking tyres
Speed control:	Twist grip speed/direction control with maximum speed governor
Instruments & Controls:	On/Off keyswitch, Power LED, Battery Gauge, Emergency Stop Switch, Horn Emergency Reverse Bar, Auto-reset Thermal cut-out, Emergency Brake Release
Capacity Std :	125 kg
Capacity HD:	225kg
Max Slope:	1 in 12



Electrodrive NZ Ltd
PO Box 58 391, Greenmount, Auckland, NZ
P: +64 9 296 3225 F: +64 9 296 3224
Web: www.electrodrive.co.nz
Email: info@electrodrive.co.nz